

# Uni Aqua

## Water Decontamination Treatment

Ready to use high-performance antimicrobial solution designed for the Continuous treatment and maintenance of dental unit waterlines (DUWLs). Powered by hydrogen peroxide and stabilized with chelators, it effectively removes and prevents biofilm, ensuring hygienic water quality in compliance with ADA and CDC recommendations.

### Active Ingredients:

- Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>) - Citric Acid - Phosphoric acid.

### Uses for Daily water decontamination & Biofilm removal:

- Use in bottle reservoir or automated waterline cleaning system
- No rinse required at low concentrations

### Storage Conditions:

- Store in original container tightly closed, away from direct sunlight
- Keep at 15–25°C
- Use within 6 months of opening

### Caution:

- For professional use only
- Do not mix with other chemicals
- Avoid contact with eyes and skin
- Keep out of reach of children

**Manufactured by:**  
UNS-EG

[www.uns-eg.com](http://www.uns-eg.com)

[info@uns-eg.com](mailto:info@uns-eg.com)

Code: 3025

Lot: AQ 001

PRO: 04-2025    EXP: 03-2027